

Cranial iD[®]

PEEK customized implant



Precise fit.

Excellent strength.

PEEK customized implant

PEEK customized implants by Stryker are designed with exacting parameters to optimize the bone-to-implant interface. Multiple fit options are available to accommodate the needs of both you and your patient. The “best in class” support and expertise of your local Stryker representative combined with the strength and fit of your implant allow you to close and reconstruct with confidence.

We have two distinct offerings based on either your need for **priority** delivery or a more **complex** surgical plan. Our comprehensive, market-leading portfolio of customized products and services allow you to tailor your surgical cases to the unique needs of your patients.

PEEK features and benefits

Precise fit

Designed and manufactured to produce a “drop-in fit”

Suture/drainage hole options

Customizable suture/drainage hole design with option for 1-6 pairs or full pattern suture hole implant; hole diameter available in either 2 or 3mm

Thickness options

Implant wall thickness provided in either 3.3, 4, 5, or 6mm to cater to patient-specific needs; tapered toward thin bone option also available

Strength



Cranial bone: 100 MPa Yield Strength¹
PEEK implant: 115 MPa Yield Strength²

Sterilization

PEEK implants are delivered non-sterile and **must** be sterilized before implantation. See IFU for sterilization requirements and re-sterilize the implant if it should ever become contaminated

Modification

If needed, the PEEK implant may be modified by burring



PEEK Priority

Priority timeline

Implants are intended to be designed, manufactured and shipped to your facility within **5 standard business days**

PEEK Complex

Standard timeline

Implants are intended to be designed, manufactured and delivered to your facility within **12 standard business days**

Host bone model

Provided in every Complex implant kit, the host bone model acts as a preoperative guide to demonstrate orientation and fit of the PEEK Complex implant

Single stage with guides

The PEEK single stage kit comes with one non-sterile implant and one non-sterile host bone. Single stage is compatible with 3D Systems guides, along with navigation

PEEK Single stage

PEEK Single stage allows the surgeon to address their patients' needs in one surgery


3D Systems' VSP® solutions

3D Systems offers surgical marking guides in cap or ring style, along with anatomical models

Navigation

Navigation can be integrated to transfer the planned resection outline to the patient's bony anatomy



Patient positioning	
Head alignment	Remain straight in neutral position.
Gantry tilt	0° gantry tilt.
 <p>No oblique angle of locator/survey lines. No gantry tilt (CT).</p>	
Scan length/Field of view (FOV)	
Scan length	For cranial defects, encompass the entire skull , including at least 2 slices superior to the skull.
FOV	For mandibular defects, encompass the entire mandible. Select FOV to include all surrounding anatomy.
Scanning process	
Patient movement	Avoid patient motion. If the scan shows motion artifacts, the scan cannot be used.
Acquisition	
Slice thickness	Maximum = 1.5 mm (1 mm preferred)
Beam collimation	Width and detector configuration necessary to achieve actual slice thickness.
Table increment	Constant table increment, no gaps. Smaller than or equal to slice thickness.
Sequential scanners	No overlap and no gap.
Single-slice helical scanners	Beam pitch = 1
Multi-slice helical scanners	Beam pitch < 1 (GE: High Quality; Toshiba: Detail)
Slice orientation	Axial slice orientation.
Algorithm (Kernel)	Bone algorithm. Warning: DO NOT post process to alter slice orientation or thickness.
Data	
Series ID	All images of a scan should be stored in one series.
File format	DICOM format. No raw data. Cone beam (CBCT) scans accepted Do not compress. Contrast not required. Inclusion of CT Viewer not recommended.
No raw data	Archive only the relevant examination in uncompressed DICOM (CD-R preferred).
Data storage	Recommendation: Save raw data for at least 14 days after scan.
Guidance for pediatric scanning	
Exposure to ionizing radiation is of particular concern in pediatric patients. Check if existing scans meets the requirements for Stryker implant design. To avoid rescanning of patients, follow the parameters given in the Stryker Scan Protocol and use reduced dose and child - sized protocols where appropriate. Stryker recommends consulting the instructions for use provided by your imaging device manufacturer, and limiting radiation dosage to the amount clinically necessary. Statutory national Diagnostic Reference Levels (DRLs) for pediatric as well as for adult CT examinations must be complied with. Limit the dose by reducing Tube Voltage (kV) and the Tube-Current-Time product (mAs), consider patient size and activate Tube Current Modulation and/or Automatic Exposure Control if applicable and indicated for pediatric patients.	

Craniomaxillofacial

References

1. J. Motherway, P. Verschueren, G. Ferre, J. Sloten, M. Gilchrist. "The Mechanical Properties of Cranial Bone: The Effect of Loading Rate and Cranial Sampling Position".
2. InVibio Limited. PEEK-Optima Polymer Typical Material Properties Brochure.
3. Liu JK, Gottfried ON, Cole CD, Dougherty WR, CouldwellWT. "MEDPOR Porous Polyethylene implant for Cranioplasty and Skull Base Reconstruction" Neurosurgery. April 2004.

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Cranial ID, iD Solutions, MEDPOR, Stryker, Stryker Orthopaedics. All other trademarks are trademarks of their respective owners or holders.

Stryker Craniomaxillofacial
 Kalamazoo, MI 49002 USA
 t: 269 389 5346
 f: 877 648 7114
 toll free: 800 962 6558

stryker.com/cmf

CCI-BR-7_Rev.None_29598
 Copyright © 2021 Stryker
 Printed in USA