

stryker®


Craniomaxillofacial



DMO

Conformable and Resorbable Collagen Membrane

DMO




Ordering Information

CDSL11	1 in x 1 in
CDSL13	1 in x 3 in
CDSL22	2 in x 2 in
CDSL33	3 in x 3 in
CDSL45	4 in x 5 in

Product Information

- Supplied sterile in double peel pouch
- Store at room temperature
- Prehydrate in sterile saline for 30 seconds
- Can be cut in the dry or hydrated state
- Apply either side down
- Ensure an overlap of 1 cm to cover the existing dura



Excellent handling characteristics

Features and Benefits

- Conformable
- Highly Purified Type I Collagen
- Resorbable
- Provides Support for Host Tissue Regeneration
- Excellent Handling Characteristics

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**. All other trademarks are trademarks of their respective owners or holders.

DuraMatrix-Onlay is manufactured by Collagen Matrix, Inc. Franklin Lakes, NJ

Part Number: 9410-400-053 Rev. A

Copyright © 2009 Stryker
Printed in USA

Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 324 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com