




MEDPOR® TSI


Designed to repair the sellar floor

MEDPOR® TSI

- MEDPOR Implant specifically designed to repair the sellar floor following a pituitary tumor resection
- Available in two sizes and configurations: original and large TSI
 - A nonporous sheet of polyethylene heat bonded to the posterior surface of the larger TSI forms a BARRIER to aid in preventing tissue ingrowth
- Easily cut with a variety of surgical instruments
- The firm nature of the material allows carving with a sharp instrument without collapsing the pore structure
- Tabs oriented on the underside of the implant to facilitate handling and placement
- Implant fixation is not required
- MEDPOR Biomaterial has a proven history of clinical success, with applications in over 400,000 cases and references in more than 350 published clinical reports in cranial, reconstructive, oculoplastic and cosmetic applications¹



CAT#	DESCRIPTION	A	B	C	D
82007	TSI	20mm	20mm	2.5mm	0.45mm
82008	TSI BARRIER - Large	40mm	40mm	2.5mm	0.73mm



1. Sales data on file at Stryker

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **MEDPOR**, **Stryker**. All other trademarks are trademarks of their respective owners or holders.

stryker®

Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 324 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com

Literature Number:
9410-400-255 Rev. None
UnDe/P.S.
Copyright © 2012 Stryker
Printed in USA