


stryker



MEDPOR TITAN[®]

3D Orbital Floor

A Stryker Exclusive Solution

MEDPOR TITAN

3D Orbital Floor

Anatomic

Implants designed using CT-scan data to approximate the anatomy of the orbital floor & medial wall to enhance the effectiveness and efficiency of reconstruction

*Packaged sterile

Proven

MEDPOR is a clinically proven material with over 400,000 procedures performed, and more than 350 published clinical reports

Smooth

MEDPOR coating minimizes sharp edges even if the plates require modification, and the superior, non-porous barrier side helps prevent tissue ingrowth along the aspect of the globe

CAT#	DESCRIPTION
81041	MEDPOR TITAN 3D OF, MTB Left Small
81042	MEDPOR TITAN 3D OF, MTB Right Small
81043	MEDPOR TITAN 3D OF, MTB Left Large
81044	MEDPOR TITAN 3D OF, MTB Right Large
01-01820	Plate Holding Forcep
01-02395	Globe Retractor, Adjustable




Plate	DESCRIPTION	A	B	C
LARGE L/R	Height of medial wall	36 mm (1.4 in.)	37 mm (1.4 in.)	17 mm (0.6 in.)
Small l/r	Height of medial wall	32 mm (1.2 in.)	35 mm (1.4 in.)	13 mm (0.5 in.)
	Profile Height (d)	1.2 mm		

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Hybrid MMF, SMARTLock, Stryker. All other trademarks are trademarks of their respective owners or holders.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 389 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com

CMF-FL-6 Rev. None
UnDe/P.S.
Copyright © 2016 Stryker
Printed in USA