



Leibinger

1.2/1.7 Bone Graft Module

1.2/1.7 Bone Graft Module

Bone Graft Module

Holds 1.2 and 1.7 diameter screws for bone grafting procedures.

Unique Features:

- SMARTLoad screw field for easy loading
- Extensive screw field to hold more screws

REF	Description
	Bone Graft Module
29-22176	Bone Graft Fixation Module
	1.2mm Self-Tapping Screws, Gold (5/Pkg)
50-12003	1.2 x 3mm
50-12004	1.2 x 4mm
50-12005	1.2 x 5mm
50-12006	1.2 x 6mm
50-12007	1.2 x 7mm
50-12008	1.2 x 8mm
50-12010	1.2 x 10mm
50-12012	1.2 x 12mm
	1.2mm Self-Drilling Screws, Silver (5/Pkg)
50-12903	1.2 x 3mm
50-12904	1.2 x 4mm
	1.4mm Emergency Screws, Purple (5/Pkg)
50-14003	1.4 x 3mm
50-14005	1.4 x 5mm
50-14007	1.4 x 7mm
50-14009	1.4 x 9mm
	Twist Drills for 1.2mm Screws
60-09504	0.9 x 46mm, 4mm working length, Stryker End
60-10506	1.0 x 50mm, 6mm working length, Stryker End
60-10512	1.0 x 54mm, 12mm working length, Stryker End
	Instruments
62-12170	Screwdriver Blade 1.2/1.7mm
62-20285	Screwdriver Handle, Universal
	1.7mm Self-Tapping Screws, Gold (5/Pkg)
50-17003	1.7 x 3mm
50-17004	1.7 x 4mm
50-17005	1.7 x 5mm
50-17006	1.7 x 6mm
50-17007	1.7 x 7mm
50-17008	1.7 x 8mm
50-17010	1.7 x 10mm
50-17012	1.7 x 12mm
	1.7mm Self-Drilling Screws, Silver (5/Pkg)
50-17903	1.7 x 3mm
50-17904	1.7 x 4mm
50-17905	1.7 x 5mm
50-17906	1.7 x 6mm
	1.9mm Emergency Screws, Purple (5/Pkg)
50-19003	1.9 x 3mm
50-19005	1.9 x 5mm
50-19007	1.9 x 7mm
50-19009	1.9 x 9mm
	Twist Drills for 1.7mm Screws
60-13505	1.35 x 50mm, 5mm working length, Stryker End
60-13506	1.35 x 50mm, 6mm working length, Stryker End
60-14008	1.4 x 79mm, 8mm working length, Stryker End
60-13512	1.4 x 54mm, 12mm working length, Stryker End
60-13570	1.3 x 70mm, 12mm working length, Stryker End

Joint Replacements

Trauma, Extremities & Deformities

Craniofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Spine

Navigation

Endoscopy

Communications

Imaging

Patient Care & Handling Equipment

EMS Equipment

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**. All other trademarks are trademarks of their respective owners or holders.

Literature Number: **9410-300-006 Rev. A**
UnDe/P.S.

Copyright © 2011 Stryker
Printed in USA

Stryker Craniofacial
Kalamazoo, MI 49002 USA
t: 269 324 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com