


AXS

Axial stability screw

**The next generation screw
for neuro fixation**




AXS

Axial stability screw


Our latest iteration of the **Universal Neuro screw** delivers on Stryker's commitment to advancing neurotechnology. A brand new head design provides surgeons with the ease-of-use they desire with the performance they demand.

- The enhanced pick-up is designed to help improve complications related to handing off a loaded screwdriver from tech to surgeon
- Self-centering feature facilitates off-axis insertion (fig.1)
- Low profile heights:
 - 4mm SD screw with blue plate - 0.55mm
 - 4mm SD screw with gold plate - 0.6mm



Cat#(set of 5) Description Dimensions(mm)

Cat#(set of 5)	Description	Dimensions(mm)
62-15035	Screwdriver blade - long	72mm
62-15036	Screwdriver blade - short	42mm
56-15034	Self-tapping screw	1.5x4mm
56-15035	Self-tapping screw	1.5x5mm
56-15036	Self-tapping screw	1.5x6mm
56-15933	Self-drilling screw	1.5x3mm
56-15934	Self-drilling screw	1.5x4mm
56-15935	Self-drilling screw	1.5x5mm
56-17334	Emergency screw	1.7x4mm



(fig.1)

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 389 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com

CMF-FL-2 Rev. 1_16004
Copyright © 2017 Stryker
Printed in USA