

AXS

Axial stability screw



1.2mm and 1.7mm

AXS

- Improved screw pick-up performance
- Improved off-axis stability during insertion (fig. 1)
- Lower profile height of screw head
- Improved insertion of self-drilling screws



(fig.1)

1.2mm Self-drilling screws

Product #	Description
56-12903	1.2x3mm
56-12904	1.2x4mm

1.2mm Self-tapping screws

Product #	Description
56-12003	1.2x3mm
56-12004	1.2x4mm
56-12005	1.2x5mm
56-12006	1.2x6mm
56-12007	1.2x7mm
56-12008	1.2x8mm
56-12010	1.2x10mm
56-12012	1.2x12mm

1.4mm Emergency screws

Product #	Description
56-14003	1.4x3mm
56-14005	1.4x5mm
56-14007	1.4x7mm
56-14009	1.4x9mm

1.2/1.7mm Screwdriver blades

Product #	Description
62-12172	42mm - short
62-12174	82mm - standard
62-12177	112mm - long

1.7mm Self-drilling screws

Product #	Description
56-17903	1.7x3mm
56-17904	1.7x4mm
56-17905	1.7x5mm
56-17906	1.7x6mm

1.7mm Self-tapping screws

Product #	Description
56-17003	1.7x3mm
56-17004	1.7x4mm
56-17005	1.7x5mm
56-17006	1.7x6mm
56-17007	1.7x7mm
56-17008	1.7x8mm
56-17010	1.7x10mm
56-17012	1.7x12mm

1.7mm Locking screws

Product #	Description
56-17504	1.7x4mm
56-17505	1.7x5mm
56-17506	1.7x6mm

1.9mm Emergency screws

Product #	Description
56-19003	1.9x3mm
56-19005	1.9x5mm
56-19007	1.9x7mm
56-19009	1.9x9mm