

# VSP® Single jaw splints

## Virtual Surgical Planning

### Surgeon and patient information

Surgeon name:			
Phone number:		Email address:	
Additional contact:			
Hospital affiliation:			

Patient name/identifier:	
Patient DOB:	
Surgery date:	
Stryker sales rep name:	

**Method for approving occlusion:** Web meeting      Case report images (must be approved prior to manufacture)

**Occlusal notes:** (ex. segmentation between 8/9, 3mm overjet, >0.5mm enameloplasty)


### Splint features:

#### Buccal labial flange

Normal (1-2.5mm flange from the buccal labial surface of the teeth)



Minimal (<1mm flange from the buccal labial surface of the teeth)



None




#### Fixation

Holes (holes placed between the maxillary teeth)



Wire groove (groove connecting holes on the underside of the splint)



None



#### Palatal support

Palatal strap (5-7mm strap spanning the width of the palate and providing support to maxillary segments)



Thick (3-5mm extending up from the occlusal surface to provide support and rigidity)



Normal (2-3mm extending up from the occlusal surface to provide support and rigidity)




None



#### Special

Palatal insert (a splint with no occlusal coverage that wraps the palate and interlocks with the final splint)



Minimal (<1mm flange from the buccal labial surface of the teeth)



### Additional notes:



Completed order form can be emailed to [vsp@3dsystems.com](mailto:vsp@3dsystems.com).

#### Manufactured by:

3D Systems  
5381 South Alkire Circle  
Littleton, Colorado 80127 USA  
t: +1 844.643.1001 (toll-free US/Canada)  
t: +1 720.643.1001  
f: +1 720.643.1009  
[www.3dsystems.com/healthcare](http://www.3dsystems.com/healthcare)  
[vsp@3dsystems.com](mailto:vsp@3dsystems.com)  
MM-649 Rev. A

#### Distributed by:

Stryker Craniomaxillofacial  
Kalamazoo, MI 49002 USA  
t: 269 389 5346  
toll free: 800 962 6558  
f: 877 648 7114

#### [stryker.com/cmfi](http://stryker.com/cmfi)

CMF-SS-5\_Rev. None\_19170  
Copyright © 2017 Stryker  
Printed in USA