

The driver made to work with
our market-leading screws

AXS™ – Axial stability screws



Enhanced pick-up
helps improve complications related to
handing off a loaded screwdriver



Low profile screw heights
4mm SD screw with blue plate - 0.55mm
4mm SD screw with gold plate - 0.6mm



Self-centering feature
facilitates off-axis insertion



Our innovative spin
on restoration

The AXS™ Neuro Driver puts you in control when placing screws during your cranial procedures. To trial the AXS™ Neuro Driver or to get more information contact your Craniomaxillofacial sales representative, or

call 800 962 6558 or visit cmf.stryker.com

Part number	Part description
PDCMF-1000-1	AXS™ Neuro Driver
PDCMF-BP-001	Battery, AXS™ Neuro Driver
PDCMF-SC-001	Storage and transport tray
JK341	1/2 size rigid container
JK388	1/2 size rigid container lid

Craniomaxillofacial

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: AXS, Stryker. All other trademarks are trademarks of their respective owners or holders.

Distributed by:
Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 389 5346
toll free: 800 962 6558
f: 877 648 7114
stryker.com/cmf

Manufactured By:
Pro-Dex, Inc.
2361 McGaw Ave.
Irvine, California 92614 USA

Copyright © 2021
Stryker Printed in USA

NF-BR-4_Rev. None_29157

stryker

stryker

AXS™ Neuro Driver

Close and restore
with **confidence**



Close and restore with **confidence**

For any case, for every patient, restoration
is in your hands.

Adaptive torque limiting functionality

The AXS Neuro Driver is programmed with a smart algorithm to sense when a screw has been inserted and stops at the optimal time. You retain full control if you need to stop sooner, simply release the button.




Continuous drill mode
Blue lights indicate AXS™ Neuro Driver is in non-torque limiting/drill mode



Auto insertion screw mode
Green lights indicate AXS™ Neuro Driver is in auto insert mode for screw insertion




Error indicator
Allows the user to troubleshoot based upon real time visual feedback



Part number	Part description
PDCMF-SC-001	Storage and transport tray, AXS™ Neuro Driver
JK341 *	½ size rigid container*
JK388	Lid ½ size rigid container

Storage and transport tray

Everything you need
in one place



Highlights

Multi-speed screw insertion for delicate procedures

Auto-insertion functionality
compatible with market-leading AXS neuro screws

Single use battery
always ready to go

Convenient reverse function for removing screws

Drill mode/screw insertion mode
for maximum versatility

