

Cranial iD[®]

PEEK Single Stage




Precise fit. Excellent **strength.**

PEEK Single Stage


Individually designed. **Personalized** care.

PEEK Single Stage allows the surgeon to address their patients' needs in one surgery. Through pre-surgical planning a predictive craniotomy is made, leading to the design and manufacturing of a patient-specific implant delivered to the hospital prior to surgery.

stryker |  3D SYSTEMS




Cap (outer marking wall)



Ring (inner marking wall)


3D Systems VSP® solutions

For PEEK Single Stage cases, 3D Systems offers surgical marking guides and anatomical models to assist with the craniotomy. Marking guides are available in either a cap or ring style.



Navigation

Navigation can be integrated into PEEK Single Stage cases. Transfer the planned resection outline to the patient's bony anatomy using compatible navigation system (STL and DICOM files).



Craniomaxillofacial

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Cranial iD, PEEK, Stryker. All other trademarks are trademarks of their respective owners or holders.

Distributed by:
Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 389 5346
toll free: 800 962 6558
f: 877 648 7114

stryker.com/cmf

Copyright © 2021 Stryker
Printed in USA

CCI-BR-6_Rev.None_29423