

stryker

VariSpeed™

Battery Powered
Screwdriver



Battery Powered Screwdriver

Enhancements

- Increased speed - maximum speed of 400rpm¹
- Increased torque - maximum torque of 1.0Nm¹
- New sensors - allow for regulation of both speed and torque
- Improved ergonomics - robust design offering enhanced grip and hand positioning

Catalog #	Description
62-50200	VariSpeed Battery Powered Screwdriver
62-50201	VariSpeed battery pack
29-15200	VariSpeed sterilization tray
62-50104	Drill Adapter - Dental
62-50105	Drill Adapter - J-Latch

Designed for drilling pilot holes and placing screws of a 1.2 mm to 2.3 mm diameter

Touch controls allow for forward and reverse motion



Battery is packaged sterile, single use

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: VariSpeed, Stryker. All other trademarks are trademarks of their respective owners or holders.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Sources: 1. VariSpeed IFU 90-02024

Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 389 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com

CMF-FL-72_Rev. None_15573
Copyright © 2017 Stryker
Printed in USA